

Catalog #: 50548

Aliases

PD1; Programmed cell death protein 1; Protein PD-1; hPD-1; CD279

Background

Gene Name: PDCD1

NCBI Gene Entry: [5133](#)

UniProt Entry: [Q15116](#)

Application Information

Molecular Weight: Predicted, 31 kDa; observed, 50,32 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD279

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD279 expression in K562 (A), U87MG (B) whole cell lysates.