

**Catalog #: 5055**

## Aliases

IKZF3; IKAROS Family Zinc Finger 3; ZNFN1A3; Aiolos; Zinc Finger Protein, Subfamily 1A, 3 (Aiolos); Zinc Finger Protein Aiolos; Zinc Finger DNA Binding Protein Aiolos; IKAROS Family Zinc Finger 3 (Aiolos); Ikaros Family Zinc Finger Protein 3; IMD84; AIO

## Background

Gene Name: IKZF3

NCBI Gene Entry: [22806](#)

UniProt Entry: [Q9UKT9](#)

## Application Information

Molecular Weight: Predicted, 58 kDa; observed, 65 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14065

Species Reactivity: Human

Applications Tested: Western blotting (WB)

## Immunogen

Recombinant protein of human IKZF3

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis using anti-IKZF3 antibody (Cat#5055). Total cell lysates (30  $\mu$ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-IKZF3 antibody (Cat#5055, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)