

Catalog #: 50598

Aliases

PTPK; Receptor-type tyrosine-protein phosphatase kappa; Protein-tyrosine phosphatase kappa; R-PTP-kappa

Background

Gene Name: PTPRK

NCBI Gene Entry: [5796](#)

UniProt Entry: [Q15262](#)

Application Information

Molecular Weight: Predicted, 162 kDa; observed, 210,162 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PTP kappa

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PTP kappa expression in HCT116 (A), A549 (B), mouse lung (C), rat brain (D) whole cell lysates.