

Catalog #: 50599

Aliases

PMP70; PXMP1; ATP-binding cassette sub-family D member 3; 70 kDa peroxisomal membrane protein; PMP70

Background

Gene Name: ABCD3

NCBI Gene Entry: [5825](#)

UniProt Entry: [P28288](#)

Application Information

Molecular Weight: Predicted, 75 kDa; observed, 75 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ABCD3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ABCD3 expression in HepG2 (A), MCF7 (B), SHSY5Y (C), SP2/0 (D), PC12 (E) whole cell lysates.