

Catalog #: 50625

Aliases

D6S218E; 40S ribosomal protein S18; Ke-3; Ke3

Background

Gene Name: RPS18

NCBI Gene Entry: [6222](#)

UniProt Entry: [P62269](#)

Application Information

Molecular Weight: Predicted, 17 kDa; observed, 18 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, pig, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RPS18

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RPS18 expression in HGC27 (A), mouse liver (B), mouse lung (C), rat liver (D), rat lung (E) whole cell lysates.