

Anti-Alpha-sarcoglycan Rabbit Polyclonal Antibody



Catalog #: 50637

Aliases

ADL; DAG2; Alpha-sarcoglycan; Alpha-SG; 50 kDa dystrophin-associated glycoprotein; 50DAG; Adhalin; Dystroglycan-2

Background

Gene Name: SGCA

NCBI Gene Entry: [6442](#)

UniProt Entry: [Q16586](#)

Application Information

Molecular Weight: Predicted, 42 kDa; observed, 43 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Alpha-sarcoglycan

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Alpha-sarcoglycan expression in HuvEc (A), Raw264.7 (B), PC12 (C) whole cell lysates.