

Catalog #: 50657

Aliases

3-oxo-5-alpha-steroid 4-dehydrogenase 1; SR type 1; Steroid 5-alpha-reductase 1; S5AR 1

Background

Gene Name: SRD5A1

NCBI Gene Entry: [6715](#)

UniProt Entry: [P18405](#)

Application Information

Molecular Weight: Predicted, 29 kDa; observed, 29 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human SRD5A1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of SRD5A1 expression in HEK293T (A), PC3 (B), A549 (C) whole cell lysates.