

Catalog #: 50673

Aliases

TADA2L; Transcriptional adapter 2-alpha; Transcriptional adapter 2-like; ADA2-like protein

Background

Gene Name: TADA2A

NCBI Gene Entry: [6871](#)

UniProt Entry: [O75478](#)

Application Information

Molecular Weight: Predicted, 51 kDa; observed, 51 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, monkey, pig, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ADA2A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ADA2A expression in HEK293T (A), A375 (B), Hela (C), mouse liver (D), mouse kidney (E), rat liver (F) whole cell lysates.