

Catalog #: 50679

Aliases

Dual specificity testis-specific protein kinase 1; Testicular protein kinase 1

Background

Gene Name: TESK1

NCBI Gene Entry: [7016](#)

UniProt Entry: [Q15569](#)

Application Information

Molecular Weight: Predicted, 67 kDa; observed, 68 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TESK1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of TSK1 expression in rat liver (A), rat muscle (B) whole cell lysates.