

Catalog #: 50682

Aliases

SML1; Trefoil factor 2; Spasmolysin; Spasmolytic polypeptide; SP

Background

Gene Name: TFF2

NCBI Gene Entry: [7032](#)

UniProt Entry: [Q03403](#)

Application Information

Molecular Weight: Predicted, 14 kDa; observed, 17 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Spasmolysin

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Spasmolysin expression in HEK293T (A), mouse liver (B), rat liver (C), rat brain (D) whole cell lysates.