

Catalog #: 50722

Aliases

CGT; UGT4; 2-hydroxyacylsphingosine 1-beta-galactosyltransferase; Ceramide UDP-galactosyltransferase; Cerebroside synthase; UDP-galactose-ceramide galactosyltransferase

Background

Gene Name: UGT8

NCBI Gene Entry: [7368](#)

UniProt Entry: [Q16880](#)

Application Information

Molecular Weight: Predicted, 61 kDa; observed, 60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human UGT8

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of UGT8 expression in mouse heart (A), mouse liver (B), rat heart (C) whole cell lysates.