

Catalog #: 50749

Aliases

ZFYVE2; Early endosome antigen 1; Endosome-associated protein p162; Zinc finger FYVE domain-containing protein 2

Background

Gene Name: EEA1

NCBI Gene Entry: [8411](#)

UniProt Entry: [Q15075](#)

Application Information

Molecular Weight: Predicted, 162 kDa; observed, 162 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human EEA1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of EEA1 expression in A431 (A), HeLa (B), THP1 (C) whole cell lysates.