

Catalog #: 50766

Aliases

CIS2; SSI2; STATI2; Suppressor of cytokine signaling 2; SOCS-2; Cytokine-inducible SH2 protein 2; CIS-2; STAT-induced STAT inhibitor 2; SSI-2

Background

Gene Name: SOCS2

NCBI Gene Entry: [8835](#)

UniProt Entry: [O14508](#)

Application Information

Molecular Weight: Predicted, 22 kDa; observed, 33 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, pig, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human SOCS2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of SOCS2 expression in HEK293T (A), Hela (B), mouse lung (C), mouse liver (D), rat liver (E) whole cell lysates.