

Catalog #: 5078

Aliases

MMAB; Metabolism Of Cobalamin Associated B; CFAP23; CblB; Methylmalonic Aciduria (Cobalamin Deficiency) CblB Type; ATP:Co(I)Rrinoid Adenosyltransferase MMAB; Cilia And Flagella Associated Protein 23; Methylmalonic Aciduria Type B Protein; ATP:Cob(I)Alamin Adenosyltransferase; Corrinoid Adenosyltransferase MMAB; Cob(I)Yrinic Acid A,C-Diamide Adenosyltransferase, Mitochondrial; Methylmalonic Aciduria (Cobalamin Deficiency) Type B; Cob(II)Yrinic Acid A,C-Diamide Adenosyltransferase; Aquocob(I)Alamin Vitamin B12s Adenosyltransferase; Cobinamide/Cobalamin Adenosyltransferase; ATP:Corrinoid Adenosyltransferase; Cob(II)Alamin Adenosyltransferase; Corrinoid Adenosyltransferase; EC 2.5.1.-; TR; Cob

Background

Gene Name: MMAB

NCBI Gene Entry: [326625](#)

UniProt Entry: [Q96EY8](#)

Application Information

Molecular Weight: Predicted, 27 kDa; observed, 23 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14205

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human MMAB

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Western Blotting (WB): 1:400-1:2,000

Note: This product is for research use only.

Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-MMAB antibody (Cat#5078). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-MMAB antibody (Cat#5078, 1:1,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM