

Catalog #: 50853

Aliases

TC21; Ras-related protein R-Ras2; Ras-like protein TC21; Teratocarcinoma oncogene

Background

Gene Name: RRAS2

NCBI Gene Entry: [22800](#)

UniProt Entry: [P62070](#)

Application Information

Molecular Weight: Predicted, 23 kDa; observed, 23 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RRAS2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RRAS2 expression in Jurkat (A), A431 (B), MCF7 (C), SP2/0 (D) whole cell lysates.