

Anti-ZC3H13 Rabbit Polyclonal Antibody



Catalog #: 50859

Aliases

KIAA0853; Zinc finger CCH domain-containing protein 13

Background

Gene Name: ZC3H13

NCBI Gene Entry: [23091](#)

UniProt Entry: [Q5T200](#)

Application Information

Molecular Weight: Predicted, 196 kDa; observed, 197 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ZC3H13

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ZC3H13 expression in Hela (A), H446 (B) whole cell lysates.