

Catalog #: 50865

Aliases

KIAA0197; NUP120; Nuclear pore complex protein Nup160; 160 kDa nucleoporin; Nucleoporin Nup160

Background

Gene Name: NUP160

NCBI Gene Entry: [23279](#)

UniProt Entry: [Q12769](#)

Application Information

Molecular Weight: Predicted, 162 kDa; observed, 160 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NUP160

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of NUP160 expression in HEK293T (A), U2OS (B), K562 (C), mouse lung (D), mouse liver (E), rat lung (F), rat liver (G) whole cell lysates.