

Catalog #: 50871

Aliases

NCKX2; Sodium/potassium/calcium exchanger 2; Na(+)/K(+)/Ca(2+)-exchange protein 2; Retinal cone Na-Ca+K exchanger; Solute carrier family 24 member 2

Background

Gene Name: SLC24A2

NCBI Gene Entry: [25769](#)

UniProt Entry: [Q9UI40](#)

Application Information

Molecular Weight: Predicted, 73 kDa; observed, 73 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NCKX2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of NCKX2 expression in HT1080 (A), MCF7 (B), Raw264.7 (C), PC12 (D) whole cell lysates.