

Catalog #: 50873

Aliases

Ras-related protein Rab-26

Background

Gene Name: RAB26

NCBI Gene Entry: [25837](#)

UniProt Entry: [Q9ULW5](#)

Application Information

Molecular Weight: Predicted, 27 kDa; observed, 28 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RAB26

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RAB26 expression in A549 (A), HepG2 (B), LO2 (C), C6 (D), mouse brain (E) whole cell lysates.