

Catalog #: 50883

Aliases

GPR104; GPR81; HCA1; Hydroxycarboxylic acid receptor 1; G-protein coupled receptor 104; G-protein coupled receptor 81

Background

Gene Name: HCAR1

NCBI Gene Entry: [27198](#)

UniProt Entry: [Q9BXC0](#)

Application Information

Molecular Weight: Predicted, 39 kDa; observed, 43 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GPR81

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GPR81 expression in HepG2 (A), HEK293T (B), BV2 (C) whole cell lysates.