

Anti-HIPK2 Rabbit Polyclonal Antibody



Catalog #: 50888

Aliases

Homeodomain-interacting protein kinase 2; hHIPk2

Background

Gene Name: HIPK2

NCBI Gene Entry: 28996

UniProt Entry: Q9H2X6

Application Information

Molecular Weight: Predicted, 130 kDa; observed, 140 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human HIPK2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HIPK2 expression in Myla2059 (A), A549 (B), SP20 (C) whole cell lysates.