

Anti-Apolipoprotein E Receptor 2 Rabbit Polyclonal Antibody



Catalog #: 50915

Aliases

APOER2; Low-density lipoprotein receptor-related protein 8; LRP-8; Apolipoprotein E receptor 2

Background

Gene Name: LRP8

NCBI Gene Entry: [7804](#)

UniProt Entry: [Q14114](#)

Application Information

Molecular Weight: Predicted, 105 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Apolipoprotein E Receptor 2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Apolipoprotein E Receptor 2 expression in HEK293T (A), H1792 (B), mouse spleen (C), rat spleen (D) whole cell lysates.