

Catalog #: 50918

Aliases

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2; NADH-ubiquinone oxidoreductase MLRQ subunit homolog; NUOMS

Background

Gene Name: NDUFA4L2

NCBI Gene Entry: [56901](#)

UniProt Entry: [Q9NRX3](#)

Application Information

Molecular Weight: Predicted, 9 kDa; observed, 10 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, bovine, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NDUFA4L2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of NDUFA4L2 expression in HepG2 (A), LO2 (B), HCT116 (C) whole cell lysates.