

Catalog #: 50922

Aliases

CD164L1; TEM1; Endosialin; Tumor endothelial marker 1; CD248

Background

Gene Name: CD248

NCBI Gene Entry: [57124](#)

UniProt Entry: [Q9HCU0](#)

Application Information

Molecular Weight: Predicted, 46,80 kDa; observed, 55,90 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD248

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD248 expression in HepG2 (A), MCF7 (B), Hela (C), H9C2 (D), mouse testis (E), rat heart (F) whole cell lysates.