

Catalog #: 50934

Aliases

SNM1B; 5' exonuclease Apollo; DNA cross-link repair 1B protein; SNM1 homolog B; SNMIB; hSNM1B

Background

Gene Name: DCLRE1B
NCBI Gene Entry: [64858](#)
UniProt Entry: [Q9H816](#)

Application Information

Molecular Weight: Predicted, 60 kDa; observed, 60,75 kDa
Clonality: Rabbit polyclonal antibody
Species Reactivity: Human, mouse
Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human DCLRE1B

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of DCLRE1B expression in A549 (A), CT26 (B) whole cell lysates.