

Catalog #: 50937

Aliases

NOGOR; Reticulon-4 receptor; Nogo receptor; NgR; Nogo-66 receptor

Background

Gene Name: RTN4R

NCBI Gene Entry: [65078](#)

UniProt Entry: [Q9BZR6](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 44-100 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RTN4R

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RTN4R expression in mouse brain (A), rat brain (B), rat spleen (C) whole cell lysates.