

Catalog #: 50945

Aliases

C20orf48; Pantothenate kinase 2 mitochondrial; hPanK2; Pantothenic acid kinase 2

Background

Gene Name: PANK2

NCBI Gene Entry: [80025](#)

UniProt Entry: [Q9BZ23](#)

Application Information

Molecular Weight: Predicted, 62 kDa; observed, 63 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PANK2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PANK2 expression in mouse muscle (A) whole cell lysates.