

Catalog #: 50946

Aliases

C15orf20; ATP-dependent DNA helicase PIF1; DNA repair and recombination helicase PIF1; PIF1/RRM3 DNA helicase-like protein

Background

Gene Name: PIF1

NCBI Gene Entry: [80119](#)

UniProt Entry: [Q9H611](#)

Application Information

Molecular Weight: Predicted, 69 kDa; observed, 70 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PIF1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PIF1 expression in HeLa (A), HEK293T (B), Raw264.7 (C), PC12 (D) whole cell lysates.