

Catalog #: 50948

Aliases

BBS14; KIAA0373; NPHP6; Centrosomal protein of 290 kDa; Cep290; Bardet-Biedl syndrome 14 protein; Cancer/testis antigen 87; CT87; Nephrocystin-6; Tumor antigen se2-2

Background

Gene Name: CEP290

NCBI Gene Entry: [80184](#)

UniProt Entry: [O15078](#)

Application Information

Molecular Weight: Predicted, 290 kDa; observed, 290 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CEP290

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CEP290 expression in HeLa (A), CCRFCM (B) whole cell lysates.