

Catalog #: 50952

Aliases

Ubiquitin carboxyl-terminal hydrolase 42; Deubiquitinating enzyme 42; Ubiquitin thioesterase 42; Ubiquitin-specific-processing protease 42

Background

Gene Name: USP42

NCBI Gene Entry: [84132](#)

UniProt Entry: [Q9H9J4](#)

Application Information

Molecular Weight: Predicted, 145 kDa; observed, 170 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human USP42

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of USP42 expression in HeLa (A), HCT116 (B), H1792 (C) whole cell lysates.