

Catalog #: 50964

Aliases

EFO1; KIAA1911; N-acetyltransferase ESCO1; CTF7 homolog 1; Establishment factor-like protein 1; EFO1p; hEFO1; Establishment of cohesion 1 homolog 1; ECO1 homolog 1; ESO1 homolog 1

Background

Gene Name: ESCO1

NCBI Gene Entry: [114799](#)

UniProt Entry: [Q5FWF5](#)

Application Information

Molecular Weight: Predicted, 94 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ESCO1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ESCO1 expression in A375 (A), SGC7901 (B), HepG2 (C), AML12 (D), PC12 (E) whole cell lysates.