

**Catalog #: 50978**

## Aliases

KIAA1983; Collagen and calcium-binding EGF domain-containing protein 1; Full of fluid protein homolog

## Background

Gene Name: CCBE1

NCBI Gene Entry: [147372](#)

UniProt Entry: [Q6UXH8](#)

## Application Information

Molecular Weight: Predicted, 44 kDa; observed, 44 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human CCBE1

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of CCBE1 expression in ACHN (A), HeLa (B), SP2/0 (C), PC12 (D) whole cell lysates.