

Catalog #: 50986

Aliases

DNA dC->dU-editing enzyme APOBEC-3A; A3A; Phorbolin-1

Background

Gene Name: APOBEC3A

NCBI Gene Entry: [100913187/200315](#)

UniProt Entry: [P31941](#)

Application Information

Molecular Weight: Predicted, 23 kDa; observed, 23 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human APOBEC3A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of APOBEC3A expression in Hela (A), HGC27 (B), mouse lung (C), mouse liver (D), rat liver (E) whole cell lysates.