

Anti-Aquaporin 12A Rabbit Polyclonal Antibody



Catalog #: 50993

Aliases

AQP12A; AQP12; AQPX2; Aquaporin-12A; AQP-12; AQP12B; Aquaporin-12B; AQP-12B

Background

Gene Name: AQP12A/AQP12B

NCBI Gene Entry: [375318/653437](#)

UniProt Entry: [Q8IXF9/A6NM10](#)

Application Information

Molecular Weight: Predicted, 31 kDa; observed, 35 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Aquaporin 12A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Aquaporin 12A expression in HEK293T (A), Hela (B), H446 (C), mouse testis (D), mouse kidney (E), rat testis (F), rat kidney (G) whole cell lysates.