

Catalog #: 51005

Aliases

CD8B1; T-cell surface glycoprotein CD8 beta chain; CD8b

Background

Gene Name: CD8B

NCBI Gene Entry: [926](#)

UniProt Entry: [P10966](#)

Application Information

Molecular Weight: Predicted, 23 kDa; observed, 45 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD8b

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD8b expression in Jurkat (A), Raw264.7 (B), H9C2 (C) whole cell lysates.