

Anti-IL-2 Rabbit Polyclonal Antibody



Catalog #: 51020

Aliases

Interleukin-2; IL-2; T-cell growth factor; TCGF; Aldesleukin

Background

Gene Name: IL2

NCBI Gene Entry: [3558](#)

UniProt Entry: [P60568](#)

Application Information

Molecular Weight: Predicted, 17 kDa; observed, 18 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human IL-2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of IL-2 expression in HEK293T (A), Hela (B), H446 (C), mouse liver (D) whole cell lysates.