

Catalog #: 51028

Aliases

ORC2L; Origin recognition complex subunit 2

Background

Gene Name: ORC2

NCBI Gene Entry: [4999](#)

UniProt Entry: [Q13416](#)

Application Information

Molecular Weight: Predicted, 65 kDa; observed, 66 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ORC2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ORC2 expression in HeLa (A), HGC27 (B), U2OS (C), A549 (D) whole cell lysates.