

Catalog #: 51029

Aliases

Phospholipid scramblase 1; PL scramblase 1; Ca(2+)-dependent phospholipid scramblase 1; Erythrocyte phospholipid scramblase; MmTRA1b

Background

Gene Name: PLSCR1

NCBI Gene Entry: [5359](#)

UniProt Entry: [O15162](#)

Application Information

Molecular Weight: Predicted, 35 kDa; observed, 38 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, pig, rabbit, sheep

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PLSCR1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PLSCR1 expression in H9C2 (A), CT26 (B), HepG2 (C), HCT116 (D), A549 (E) whole cell lysates.