

Catalog #: 51032

Aliases

PHOG; Short stature homeobox protein; Pseudoautosomal homeobox-containing osteogenic protein; Short stature homeobox-containing protein

Background

Gene Name: SHOX

NCBI Gene Entry: [6473](#)

UniProt Entry: [O15266](#)

Application Information

Molecular Weight: Predicted, 32 kDa; observed, 32 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken, dog

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human SHOX

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of SHOX expression in mouse muscle (A), rat muscle (B) whole cell lysates.