

Anti-Protachykinin 1 Rabbit Polyclonal Antibody



Catalog #: 51037

Aliases

NKA; NKNA; TAC2; Protachykinin-1; PPT

Background

Gene Name: TAC1

NCBI Gene Entry: [6863](#)

UniProt Entry: [P20366](#)

Application Information

Molecular Weight: Predicted, 15 kDa; observed, 15 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Protachykinin 1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Protachykinin 1 expression in mouse spleen (A), rat spleen (B) whole cell lysates.