

Catalog #: 51043

Aliases

HSS; KIAA0516; MAPK8IP4; SYD1; C-Jun-amino-terminal kinase-interacting protein 4; JIP-4; JNK-interacting protein 4; Cancer/testis antigen 89; CT89; Human lung cancer oncogene 6 protein; HLC-6; JNK-associated leucine-zipper protein; JLP; Mitogen-activated protein kinase 8-interacting protein 4; Proliferation-inducing protein 6; Protein highly expressed in testis; PHET; Sperm surface protein; Sperm-associated antigen 9; Sperm-specific protein; Sunday driver 1

Background

Gene Name: SPAG9

NCBI Gene Entry: [9043](#)

UniProt Entry: [O60271](#)

Application Information

Molecular Weight: Predicted, 146 kDa; observed, 190 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human JIP4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of JIP4 expression in AML12 (A), U87MG (B), EC9706 (C), SGC7901 (D) whole cell lysates.