

Catalog #: 51046

Aliases

YY1-associated factor 2

Background

Gene Name: YAF2

NCBI Gene Entry: [10138](#)

UniProt Entry: [Q8IY57](#)

Application Information

Molecular Weight: Predicted, 19 kDa; observed, 19 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human YAF2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of YAF2 expression in Jurkat (A), Raw264.7 (B), H9C2 (C) whole cell lysates.