

Catalog #: 51062

Aliases

BET; Delta and Notch-like epidermal growth factor-related receptor

Background

Gene Name: DNER

NCBI Gene Entry: [92737](#)

UniProt Entry: [Q8NFT8](#)

Application Information

Molecular Weight: Predicted, 78 kDa; observed, 110 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human BET

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of BET expression in rat kidney (A), SP20 (B), C6 (C), SHSY5Y (D) whole cell lysates.