

Anti-Collagen 26 alpha 1 Rabbit Polyclonal Antibody



Catalog #: 51064

Aliases

EMID2; EMU2; Collagen alpha-1(XXVI) chain; Alpha-1 type XXVI collagen; EMI domain-containing protein 2; Emilin and multimerin domain-containing protein 2; Emu2

Background

Gene Name: COL26A1

NCBI Gene Entry: [136227](#)

UniProt Entry: [Q96A83](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 43 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Collagen 26 alpha 1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Collagen 26 alpha 1 expression in HEK293T (A), A549 (B), mouse kidney (C), rat kidney (D), rat liver (E) whole cell lysates.