

Catalog #: 51065

Aliases

RGM domain family member B; DRG11-responsive axonal guidance and outgrowth of neurite; DRAGON

Background

Gene Name: RGMB

NCBI Gene Entry: [285704](#)

UniProt Entry: [Q6NW40](#)

Application Information

Molecular Weight: Predicted, 47 kDa; observed, 52 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RGMB

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RGMB expression in mouse brain (A), rat brain (B), HeLa (C), HCT116 (D), A2780 (E), PC3 (F), SGC7901 (G) whole cell lysates.