

Anti-Beta-3 Adrenergic Receptor Rabbit Polyclonal Antibody



Catalog #: 51075

Aliases

ADRB3R; B3AR; Beta-3 adrenergic receptor; Beta-3 adrenoreceptor; Beta-3 adrenoceptor

Background

Gene Name: ADRB3

NCBI Gene Entry: [155](#)

UniProt Entry: [P13945](#)

Application Information

Molecular Weight: Predicted, 43 kDa; observed, 44,68 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Beta-3 Adrenergic Receptor

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Beta-3 Adrenergic Receptor Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Beta-3 Adrenergic Receptor expression in HEK293T (A), Hela (B), HepG2 (C), A2780 (D), mouse kidney (E), rat liver (F) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM