

Catalog #: 51083

Aliases

Sodium/potassium-transporting ATPase subunit beta-3; Sodium/potassium-dependent ATPase subunit beta-3; ATPB-3; CD298

Background

Gene Name: ATP1B3

NCBI Gene Entry: [483](#)

UniProt Entry: [P54709](#)

Application Information

Molecular Weight: Predicted, 31 kDa; observed, 40 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD298

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD298 expression in HEK293T (A), Hela (B), A2780 (C), mouse heart (D), mouse liver (E), rat heart (F), rat liver (G) whole cell lysates.