

## Catalog #: 5118

### Aliases

IRAG2; Inositol 1,4,5-Triphosphate Receptor Associated 2; JAW1; LRMP; Lymphoid-Restricted Membrane Protein; Lymphoid Restricted Membrane Protein; Protein Jaw1

### Background

Gene Name: IRAG2

NCBI Gene Entry: [4033](#)

UniProt Entry: [Q12912](#)

### Application Information

Molecular Weight: Predicted, 16 kDa; observed, 16 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15801

Species Reactivity: Human

Applications Tested: Western blotting (WB)

### Immunogen

Recombinant protein of human IRAG2

### Isotype

Mouse IgG1

### Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### Storage

Store at -20 °C for one year.

### Recommended Dilutions

Western Blotting (WB): 1:400-1:2,000

**Note:** This product is for research use only.

### Validation Data

---


#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis using anti-IRAG2 antibody (Cat#5118). Total cell lysates (30  $\mu$ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-IRAG2 antibody (Cat#5118, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)