

Anti-Apolipoprotein B Receptor Rabbit Polyclonal Antibody



Catalog #: 51219

Aliases

APOB48R; Apolipoprotein B receptor; Apolipoprotein B-100 receptor; Apolipoprotein B-48 receptor; Apolipoprotein B48 receptor; apoB-48R

Background

Gene Name: APOBR

NCBI Gene Entry: [55911](#)

UniProt Entry: [Q0VD83](#)

Application Information

Molecular Weight: Predicted, 115 kDa; observed, 100 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Apolipoprotein B Receptor

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Apolipoprotein B Receptor expression in SP20 (A), PC12 (B), SHSY5Y (C), MCF7 (D), H460 (E) whole cell lysates.