

Catalog #: 51259

Aliases

Coagulation factor XI; FXI; Plasma thromboplastin antecedent; PTA

Background

Gene Name: F11

NCBI Gene Entry: [2160](#)

UniProt Entry: [P03951](#)

Application Information

Molecular Weight: Predicted, 70 kDa; observed, 72 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Factor XI LC

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Factor XI LC expression in Jurkat (A) whole cell lysates.